

Features:

DFS-Taq "Hot" is the "Hot-Start PCR" version of DFS Taq DNA Polymerase that provides a new formula in buffers and additives to prevent failures in PCR-applications where inhibitors (e.g. proteins, fat or PS) reduce the performance.

The Taq Polymerase is blocked with monoclonal Antibodies. Both components, the **Taq Polymerase and the Antibody is available, separately** at GeneON's shop.

The robust enzyme is well suited for sensitive experiments using random primers or bacterial templates. Because of the high sensitivity less than 6 molecules can be detected.

- reliable and reproducible quantification in qPCR
- perfect for real time PCR
- especially for diagnostic purposes
- reaction set-up at room temperature
- activation of enzyme during first heating
- no change or optimization of protocol necessary
- high specificity, reduced primer mismatch or dimers

Application:

Instead of conventionally purified Taq-DNA Polymerase for sensitive PCR reactions, for the detection of bacterial DNA or for applications where inhibitors decrease the performance of regular polymerases.

- Hot start PCR
- Real time PCR
- Amplification of complex genomic and cDNA templates
- Multiplex PCR
- High specificity PCR

Reaction conditions:

Same as for conventionally purified Taq-DNA Polymerase.

Concentration

5 units/μl supplied in 10 mM KPO₄ (pH 7.4 at 25°C), 0.1 mM EDTA, 0.1% Tween 20, 0.1% Triton-X 100 and 50 % (v/v) glycerol.

Unit Definition

One unit is defined as the amount of enzyme required to incorporate 10 nmoles of dNTP into acid-insoluble material in 30 min at 74°C.

Activity assay: 50 mM Tris-HCl (pH 8.0 at 25°C), 50 mM NaCl, 10 mM MgCl₂, 200 µM dATP, 200 µM dCTP, 200 µM dGTP, 50 µM [3H] dTTP, 0,25 mg/ml activated calf thymus DNA.

Storage conditions: -20°C in 50 mM Tris-HCl (pH 8.0 at 25°C), 50 mM NaCl, 0.1 mM EDTA, 1 mM DTT, 50% glycerol and 1% triton X-100.

Quality control: Endo-, exodeoxyribonucleases, ribonucleases free. Free of bacterial DNA-traces

Components	Volume per reaction
10X reaction buffer	5 µl
100 mM MgCl ₂	Optional
dNTP-Mix (40mM)	1.0 µl
Up-stream primer (10 µM stock)	0.5 - 2.5 µl
Down-stream primer (10 µM stock)	0.5 - 2.5 µl
Template DANN	0.1 - 15 ng/ml plasmid DANN 1-10 µg/ml genomic DANN
DFS-Hot Taq DANN (5 u/µl)	0.2 - 1.0 µl
Sterile dest. Water (molecular grade)	Up to 50 µl total reaction volume

Note:

- vortex all solutions carefully before using
- add the enzyme after Template DNA
- may you have to optimize the MgCl₂ concentration for best result

General Thermo-Cycler protocol:

Step	Time	Temperature
Initial denaturation	2-5min	94-95°C
25-30 Cycles:		
Denaturation	10-25 sec	94-95°C
Annealing	10-25 sec	55-65°C
Extension	60 sec	72°C per 1kb
Final extension	5 min	72°C

Note: In case of low amount of DNA template, additionally cycles may be used

Order information

Prod. No.	Description	Quantity
N9150	DFS Taq "HOT" PLUS DNA Polymeras	500 Units
N9152	DFS Taq "HOT" PLUS DNA Polymeras	5 x 500 Units
N9154	DFS Taq "HOT" PLUS DNA Polymeras	20x 500 Units