

## dNTP Set 4x 0.2ml (4x100 mM)

Prod. No.: 1910011

### Applications of dNTP-mix and dNTP sets:

- all molecular biology applications
- PCR/qPCR/real-time PCR
- Reverse transcription
- DNA labelling
- DNA sequencing

### Description of dNTP-mix and dNTP set:


dNTP-sets from bioswisstec have highest PCR Grade and contains four single tubes of dATP, dCTP, dGTP and dTTP supplied as aqueous solutions at pH 8.5. The dNTP-mixes contains an ideal mixture of dNTP's with 10 mM of each nucleotide.

### Quality control for dNTP mixes and dNTP set:

- free of human bacterial DNA
- 18 kb PCR amplification with 20 pg of template or less
- RT-PCR: 600 bp fragment with 20 pg template or less
- no DNase, RNase or Nicking activities
- tested for Inorganic species
- Swiss technology at ISO 9001:2000 certificated company under licence

### Components of dNTP-mixes and dNTP-sets:

<b>dATP</b> 2'-Deoxyadenosine 5'-triphosphate, sodium salt CAS Number: 1927-31-7 Formula: $C_{10}H_{13}N_5O_{12}P_3$ (Anion) Molecular weight: $488.16 \text{ g}\cdot\text{mol}^{-1}$ Density: 100 mM	<b>dCTP</b> 2'-Deoxycytidine 5'-triphosphate, sodium salt CAS Number: 102783-51-7 Formula: $C_9H_{13}N_3O_{13}P_3$ (Anion) Molecular weight: $464.13 \text{ g}\cdot\text{mol}^{-1}$ Density: 100 mM
<b>dGTP</b> 2'-Deoxyguanosine 5'-triphosphate, sodium salt CAS Number: 93919-41-6 Formula: $C_{10}H_{13}N_5O_{13}P_3$ (Anion) Molecular weight: $504.16 \text{ g}\cdot\text{mol}^{-1}$ Density: 100 mM	<b>dTTP</b> 2'-Deoxythymidine 5'-triphosphate, sodium salt CAS Number: 18423-43-3 Molecular formula: $C_{10}H_{14}N_2O_{14}P_3$ (Anion) Molecular weight: $479.14 \text{ g}\cdot\text{mol}^{-1}$ Density: 100 mM



RP-HPLC purity test  
Date: January 2010

**Storage:** store at - 20°C up to 24 months or short term storage at room temperature

**Shipping:** with blue ice

Prod.-No	Description	Quantity
1910011	set of 4 dNTP's	4 x 0.2 ml (100 mM each)
1910012	set of 4 dNTP's	4 x 1.0 ml (100 mM each)

## dNTP Set 4x 1ml (4x100 mM) Prod. No.: 1910012

### Applications of dNTP-mix and dNTP sets:

- all molecular biology applications
- PCR/qPCR/real-time PCR
- Reverse transcription
- DNA labelling
- DNA sequencing

### Description of dNTP-mix and dNTP set:


dNTP-sets from bioswisstec have highest PCR Grade and contains four single tubes of dATP, dCTP, dGTP and dTTP supplied as aqueous solutions at pH 8.5. The dNTP-mixes contains an ideal mixture of dNTP's with 10 mM of each nucleotide.

### Quality control for dNTP mixes and dNTP set:

- free of human bacterial DNA
- 18 kb PCR amplification with 20 pg of template or less
- RT-PCR: 600 bp fragment with 20 pg template or less
- no DNase, RNase or Nicking activities
- tested for Inorganic species
- Swiss technology at ISO 9001:2000 certificated company under licence

### Components of dNTP-mixes and dNTP-sets:

<b>dATP</b> 2'-Deoxyadenosine 5'-triphosphate, sodium salt CAS Number: 1927-31-7 Formula: $C_{10}H_{13}N_5O_{12}P_3$ (Anion) Molecular weight: $488.16 \text{ g}\cdot\text{mol}^{-1}$ Density: 100 mM	<b>dCTP</b> 2'-Deoxycytidine 5'-triphosphate, sodium salt CAS Number: 102783-51-7 Formula: $C_9H_{13}N_3O_{13}P_3$ (Anion) Molecular weight: $464.13 \text{ g}\cdot\text{mol}^{-1}$ Density: 100 mM
<b>dGTP</b> 2'-Deoxyguanosine 5'-triphosphate, sodium salt CAS Number: 93919-41-6 Formula: $C_{10}H_{13}N_5O_{13}P_3$ (Anion) Molecular weight: $504.16 \text{ g}\cdot\text{mol}^{-1}$ Density: 100 mM	<b>dTTP</b> 2'-Deoxythymidine 5'-triphosphate, sodium salt CAS Number: 18423-43-3 Molecular formula: $C_{10}H_{14}N_2O_{14}P_3$ (Anion) Molecular weight: $479.14 \text{ g}\cdot\text{mol}^{-1}$ Density: 100 mM



RP-HPLC purity test  
Date: January 2010

**Storage:** store at - 20°C up to 24 months or short term storage at room temperature

**Shipping:** with blue ice

Prod.-No	Description	Quantity
1910011	set of 4 dNTP's	4 x 0.2 ml (100 mM each)
1910012	set of 4 dNTP's	4 x 1.0 ml (100 mM each)