

## Boeco Disposable Cuvettes



**BRA 759007**

**Disposable Macro Cuvettes, PS**

Filling volume: min 2,5 ml / max. 4,5 ml from 10.000 pcs.  
grouped by mold cavity number, neutral packing  
Wavelength: from 340 to 900 nm  
packing: 100 pcs/box, 1000 pcs./carton

**60.- /1000 pcs.**

**53.- /1000 pcs.**



**BRA 759017**

**Disposable Semi-Micro Cuvettes, PS**

Filling volume: min 1,5 ml / max. 3,0 ml from 10.000 pcs.  
grouped by mold cavity number, neutral packing  
Wavelength: from 340 to 900 nm  
packing: 100 pcs/box, 1000 pcs./carton

**60.- /1000 pcs.**

**53.- /1000 pcs.**

**BRA 759170**

**BRAND UV Macro cuvette**

Filling volume: min 2,5 ml / max. 4,5 ml from 5.000 pcs.  
grouped by mold cavity number  
Wavelength: from 220 to 900 nm  
packing: pack of 100 pcs.

**35.-- /100 pcs.**

**29.- /100 pcs.**

**BRA 759150**

**BRAND UV Semi-Micro cuvette**

Filling volume: min 1,5 ml / max. 3,0 ml from 5.000 pcs.  
grouped by mold cavity number  
Wavelength: from 220 to 900 nm  
packing: pack of 100 pcs.

**35.-- /100 pcs.**

**29.-- /100 pcs.**



**BRA 759220**

**BRAND UV Micro cuvette, center height 15 mm**

Filling volume: min 70 µl / max. 550 µl from 4.000 pcs.  
grouped by mold cavity number  
for determination of proteins, DNA, RNA  
Wavelength: from 220 to 900 nm  
packing: pack of 100 pcs.

**22.- /100 pcs.**

**19.50 /100 pcs.**

**Minimum 1000 pcs.**

**Minimum Order disposable cuvettes 300,00**

## Optical glass and Quartz Cells

### BOECO Cells



<b>BOE 104010</b>	<b>BOECO Macro Cell, with PTFE lid</b> Light path 10 mm. <b>Optical glass</b> 360 nm - 2500 nm 2 pcs. / pack	from 2 pcs.	<b>23.-- / pc.</b>
		from 20 pcs.	<b>21.-- / pc.</b>
		from 50 pcs.	<b>19.-- / pc.</b>
<b>BOE 204010</b>	<b>BOECO Macro Cell, with PTFE lid</b> Light path 10 mm. <b>Quartz glass</b> 190 nm - 2500 nm 2 pcs. / pack	from 2 pcs.	<b>39.-- / pc.</b>
		from 20 pcs.	<b>34.-- / pc.</b>
		from 50 pcs.	<b>32.-- / pc.</b>

### Hellma Cells



<b>HEL 6030-10-10</b>	<b>Macro Cell 6030-OG, 45 x 12,5 mm, without lid</b> Light path 10 mm. Optical glass 360 nm - 2500 nm		<b>29.-- / pc.</b>
<b>HEL 6030-UV-10-531</b>	<b>Macro Cell 6030-UV, 45 x 12,5 mm, without lid</b> Light path 10 mm . Quartz HERASIL 230 nm - 2500 nm		<b>76.-- / pc.</b>
<b>HEL 100-10-20</b>	<b>Macro Cell, 100-OS 45 x 12,5 mm, with PTFE lid</b> Light path 10 mm . Special optical glass 320 nm - 2500 nm		<b>36.-- / pc.</b>
<b>HEL 100-20-20</b>	100-OS with 20 mm light path		<b>66.-- / pc.</b>
<b>HEL 100-50-20</b>	100-OS with 50 mm light path		<b>88.- / pc.</b>
<b>HEL 100-100-20</b>	100-OS with 100 mm light path		<b>163.-- / pc.</b>
<b>HEL 100-10-40</b>	<b>Macro Cell 100-QS, 45 x 12,5 mm, with PTFE lid</b> Light path 10 mm. Quartz SUPRASIL 200 nm - 2500 nm		<b>78.-- / pc.</b>
<b>HEL 100-20-40</b>	100-QS with 20 mm light path		<b>99.-- / pc.</b>
<b>HEL 100-50-40</b>	100-QS with 50 mm light path		<b>135.- / pc.</b>
<b>HEL 100-100-40</b>	100-QS with 100 mm light path		<b>219.-- / pc.</b>
<b>HEL 6040-10-10</b>	<b>Semi-Micro Cell 6040-OG, 45 x 12,5 mm, without lid</b> Light path 10 mm . Optical glass 360 nm - 2500 nm		<b>49.-- / pc.</b>
<b>HEL 6040-UV-10-531</b>	<b>Semi-Micro Cell 6040-UV, 45 x 12,5 mm, without lid</b> Light path 10 mm. Quartz HERASIL 230 nm - 2500 nm		<b>89.-- / pc.</b>
<b>HEL 105-202-15-40</b>	<b>Ultra Micro Cell 105-202-QS, 45 x 12,5 mm, with PTFE lid</b> Quartz SUPRASIL 200 nm - 2500 nm Light path 10 mm. Centre height 15 mm Chamber volume 50 µl, Filling Volume 70 µl		<b>219.-- / pc.</b>

**Minimum Order disposable cuvettes 300,00 all prices in CHF**  
**All prices are to be understood, net ex factory Schaffhausen**