

Centrifuge Tubes Quality Falcon Typ

Centrifuge Tubes by Bioswisstec® incl. 1 free Rack as stocks last

- ✓ Pyrogen free by LAL Test
- ✓ Free of detectable DNase, RNase and human DNA
- ✓ Free of detectable endotoxines, free of cytotoxic substances
- ✓ Raw material USP 6 (medical grade)
- ✓ Highest endurance in all applications
- ✓ Sterile, Gamma sterilised SAL 10⁻⁶
- ✓ Zip bag closure easy and secure reclosing
- ✓ Screw cap, aerosole tight in refreshing Lime green made out of HDPE (High density polyethylene)
- ✓ Tube, made out of polypropylene with high thermic and chemical resistance, for storage of chemical and biological samples
- ✓ Printing precision +/- 1.5%
- ✓ Manufactured with highly precise, diamond polished Swiss injection molds
- ✓ Lot number system: traceability guaranteed

Learn more



53015 Centrifuge Tubes 15 ml

| Specification | Dimensions | Packaging Unit | Price / Box | Price / >4 boxes |
|--|------------|----------------|-------------|------------------|
| PP, with HDPE cap, conical, graduated, with white labeling field | 17/120 mm | 50 pcs/bag | 56.50 | 52.— |
| | | 500 pcs/box | | |

53050 Centrifuge Tubes 50 ml

| Specification | Dimensions | Packaging Unit | Price / Box | Price / >4 boxes |
|--|------------|----------------|-------------|------------------|
| PP, with HDPE cap, conical, graduated, with white labeling field | 30/115mm | 25 pcs/bag | 58.50 | 56.— |
| | | 500 pcs/box | | |

53250 Centrifuge Tubes 50 ml skirted

| Specification | Dimensions | Packaging Unit | Price / Box | Price / >4 boxes |
|---|------------|----------------|-------------|------------------|
| PP, with HDPE cap, conical, graduated, with white marking field, selfstanding (skirted) | 30/115mm | 25 pcs/bag | 62.80 | 59.— |
| | | 500 pcs/box | | |

Price in CHF, ex works Schaffhausen, excl. 8% VAT, free delivery from CHF 400.—