


State of matter	Liquid
additional product description	filtered, autoclaved, sterile
EC-No.:	205-358-3
Storage:	RT
LGK:	10 - 13
Hazard pictogram(s)	
Hazard statement(s)	H373
Precautionary statement(s)	P260-P314-P501
Signal word	Warning
WGK:	2

Specification

DNases/RNases/Proteases	not detectable
pH (20°C)	8.0 ± 0.1
Composition:	
EDTA · Na ₂ · 2H ₂ O	186.12 g/L
NaOH (10 M)	approx. 50 ml/L

Order Information

Prod. No.	Description	Quantity
CE9141	EDTA solution pH 8.0 (0.5 M)	500 ml