

<b>Synonym</b>	<b>β-Nicotinamide Adenine Dinucleotide</b>
State of matter	Solid
Formula	C <sub>21</sub> H <sub>27</sub> N <sub>7</sub> O <sub>14</sub> P <sub>2</sub>
M	663.43 g/mol
CAS-No.:	53-84-9
EC-No.:	200-184-4
Storage:	-20°C
LGK:	10 - 13
WGK:	1
	electrostatic
<b>Specification</b>	
Assay (enzym.)	min. 97 %
Water (K.F.)	max. 10.0 %
UV Absorption:	
A250/A260 (pH 7.5)	0.83 ± 0.02
A280/A260 (pH 7.5)	0.20 ± 0.02

## Order Information

<b>Prod. No.</b>	<b>Description</b>	<b>Quantity</b>
CE9196	NAD	1000 mg
CE9197	NAD	5000 mg