

Synonym	$\beta$ -Nicotinamide Adenine Dinucleotide
State of matter	Solid
Formula	$C_{21}H_{27}N_7Na_2O_{14}P_2$
M	709.41 g/mol
CAS-No.:	606-68-8
EC-No.:	210-123-3
Storage:	-20°C
LGK:	10 - 13
WGK:	1
	electrostatic
Specification	
Assay (enzym.)	min. 95 %
Water (K.F.)	max. 8.0 %
UV Absorption:	
A250/A260 (pH 10)	$0.82 \pm 0.03$
A280/A260 (pH 10)	$0.23 \pm 0.02$

## Order Information

Prod. No.	Description	Quantity
CE9198	NADH - Disodium salt	1 g
CE9199	NADH - Disodium salt	5 g