

State of matter	Solid
usage	Water-soluble part of autolysed yeast. Rich of ammonia and vitamins, especially of water-soluble vitamin B-complexes; suitable for addition to media.
EC-No.:	232-387-9
Storage:	RT
WGK:	1
Specification:	
DNases/RNases/Proteases	not detectable
Ash	max. 14 %
Total N	11 - 12 %
pH (1 %; H ₂ O)	6.8 - 7.2 (25°C)
Heavy metals (as Pb)	max. 0.0001 %
Loss on drying	max. 5 %
Coagulated protein	not detectable
Sodium chloride	max. 0.3 %

Order Information

Prod. No.	Description	Quantity
CE9169	Yeast extract powder	500 g