

Specification

Product #: L 1803

Description: Trypsin/EDTA, 100 ml

Specification: Trypsin 0,05 %/EDTA 0,02 % in PBS, w/o: Ca²⁺ and Mg²⁺

Storage: -20°C

Testings

Acceptance Criteria

pH Value

Osmolality

Incubation at 32°C

Incubation at 20°C

Mycoplasma

Porcine parvovirus PPV

Performance

Accepted Result

7.3 - 7.5

260 - 290 mOsm/kg

Sterile

Sterile

Negative

Negative

Cell detachment

Composition

Inorganic Salts

Di-Sodium hydrogen phosphate	1,150.00 mg/L
Potassium chloride	200.00 mg/L
Potassium dihydrogen phosphate	200.00 mg/L
Sodium chloride	8,000.00 mg/L

Other Components

EDTA-Na ₂	200.00 mg/L
Trypsin	500.00 mg/L