

Centrifuge Tubes Quality Falcon Typ

Centrifuge Tubes by Bioswisstec® now even better with Zipbag and sealed bags

- ✓ RCF tested, 53015: 15000 fixed bucked rotor with Methanol
53050: 15000 fixed bucked rotor with Methanol
53250: 3500 fixed bucked rotor with Methanol
https://en.wikipedia.org/wiki/Centrifugal_force
- ✓ Pyrogen free by LAL Test
- ✓ Free of detectable DNase, RNase and human DNA
- ✓ Free of detectable endotoxines, free of cytotoxic substances
- ✓ Raw material USP 6 (medical grade)
- ✓ Highest endurance in all applications
- ✓ Sterile, Gamma sterilised SAL 10⁻⁶
- ✓ Zip bag closure easy and secure reclosing
- ✓ Screw cap, aerosole tight in refreshing Lime green made out of HDPE (High density polyethylene)
- ✓ Tube, made out of polypropylene with high thermic and chemical resistance, for storage of chemical and biological samples
- ✓ Printing precision +/- 1.5%
- ✓ Manufactured with highly precise, diamond polished Swiss injection molds
- ✓ Lot number system: traceability guaranteed

Learn more



53015 Centrifuge Tubes 15 ml

Specification	Dimensions	Packaging Unit	Price / Box	Price / >4 boxes
PP, with HDPE cap, conical, graduated, with white labeling field	17/120 mm	50 pcs/bag 500 pcs/box	56.50	52.—

53050 Centrifuge Tubes 50 ml

Specification	Dimensions	Packaging Unit	Price / Box	Price / >4 boxes
PP, with HDPE cap, conical, graduated, with white labeling field	30/115mm	25 pcs/bag 500 pcs/box	58.50	56.—

53250 Centrifuge Tubes 50 ml skirted

Specification	Dimensions	Packaging Unit	Price / Box	Price / >4 boxes
PP, with HDPE cap, conical, graduated, with white marking field, selfstanding (skirted)	30/115mm	25 pcs/bag 500 pcs/box	62.80	59.—

Price in CHF, ex works Schaffhausen, excl. 7.7% VAT, free delivery from CHF 400.—