

## Specification

**Product #:** A 1201

**Description:** Demecolcin solution 10 µg/ml, 10 ml

**Specification:** Demecolcin solution 10 µg/ml

**Storage:** -20°C

## Testings

### Acceptance Criteria

pH Value

Osmolality

Incubation at 32°C

Incubation at 20°C

### Accepted Result

as reported

as reported

Sterile

Sterile

## Composition

### Inorganic Salts

di-Sodium hydrogen phosphate anhydrous	1,150.0000 mg/L
Potassium chloride	200.0000 mg/L
Potassium dihydrogen phosphate	200.0000 mg/L
Sodium chloride	8,000.0000 mg/L

### Other Components

Demecolcine	10.0000 mg/L
-------------	--------------