

Specification

Product #: K 1501

Description: EDTA 1 % in PBS, w/o: Ca²⁺ and Mg²⁺ , 100 ml

Specification: EDTA 1 % in PBS, w/o: Ca²⁺ and Mg²⁺

Storage: 2°C - 8°C

Testings

Acceptance Criteria

pH Value

Osmolality

Incubation at 32°C

Incubation at 20°C

Accepted Result

7.0 - 7.3

n.a. mOsm/kg

Sterile

Sterile

Composition

Inorganic Salts

Di-Sodium hydrogen phosphate	1,150.0000 mg/L
Potassium chloride	200.0000 mg/L
Potassium dihydrogen phosphate	200.0000 mg/L
Sodium chloride	8,000.0000 mg/L

Other Components

EDTA	10,000.0000 mg/L
------	------------------