

Specification

Product #: L 1411

Description: SEcoll human, density/Dichte 1,077 g/ml, 100 ml

Specification: density/Dichte 1,077 g/ml

Storage: 2° - Room Temperature

Testings

Acceptance Criteria

Appearance

pH Value

Osmolality

Density

Accepted Result

Clear, colourless

6.6 - 7.5

255 - 275 mOsm/kg

1.076 - 1.078 g/ml

Composition

Diatrizoic Acid

86,900.0000 mg/L

Polysucrose

57,000.0000 mg/L