

## Specification

**Product #:** L 1412

**Description:** SEcoll human, 500 ml

**Specification:** density/Dichte 1,077 g/ml

**Storage:** 2° - Room Temperature

## Testings

### Acceptance Criteria

Appearance

pH Value

Osmolality

Density

### Accepted Result

Clear, colourless

6.6 - 7.5

255 - 275 mOsm/kg

1.076 - 1.078 g/ml

## Composition

Diatrizoic Acid

86,900.0000 mg/L

Polysuccrose

57,000.0000 mg/L