

Specification

Product #: L 1415

Description: Biocoll, 500ml

Specification: Biocoll-Separating Solution, desity 1.10 g/ml isoton

Storage: 2°C - 8°C

Testings

Acceptance Criteria

Appearance

pH Value

Osmolality

Density

Endotoxin

Incubation at 32°C

Incubation at 20°C

Accepted Result

Clear, colourless

6.5 - 7.5

290 - 310 mOsm/kg

1.076 - 1.077 g/ml

as reported

Sterile

Sterile

Composition

Diatrizoic Acid	86,900.0000 mg/L
-----------------	------------------

Polysuccrose	57,000.0000 mg/L
--------------	------------------