

Specification

Product #: L 1824

Description: Trypsin 0,25 %/EDTA 0,02 % in PBS, w/o: Ca²⁺ and Mg²⁺ , 100 ml

Specification: Trypsin 0,25 %/EDTA 0,02 % in PBS, w/o: Ca²⁺ and Mg²⁺

Storage: -20°C

Testings

Acceptance Criteria

pH Value

Osmolality

Incubation at 32°C

Incubation at 20°C

Accepted Result

7.3 - 7.5

as reported

Sterile

Sterile

Composition

Inorganic Salts

Di-Sodium hydrogen phosphate	1,150.0000 mg/L
Potassium chloride	200.0000 mg/L
Potassium dihydrogen phosphate	200.0000 mg/L
Sodium chloride	8,000.0000 mg/L

Other Components

EDTA-Na ₂	200.0000 mg/L
Trypsin	2,500.0000 mg/L