

Specification

Product #: L 1841

Description: Trypsin (10x) 0,05 %/EDTA 4 Na 5.3 mM w/o Phenol red w/o Ca and Mg, 100 ml

Specification: Trypsin (10x) 0,05 %/EDTA 4 Na 5.3 mM w/o Phenol red w/o Ca and Mg

Storage: -20°C

Testings

Acceptance Criteria

pH Value

Osmolality

Incubation at 32°C

Incubation at 20°C

Performance

Accepted Result

7.3 - 7.5

as reported

Sterile

Sterile

Cell detachment

Composition

Inorganic Salts

Sodium chloride

8,500.0000 mg/L

Other Components

EDTA-Na₄ x 2H₂O

2,000.0000 mg/L

Trypsin

5,000.0000 mg/L