

Specification

Product #: L 1843

Description: Trypsin (10x) 0,5%/EDTA 0,2% in PBS, w/o: Ca²⁺ and Mg²⁺, 100 ml

Specification: Trypsin (10x) 0,5%/EDTA 0,2% in PBS, w/o: Ca²⁺ and Mg²⁺

Storage: -20°C

Testings

Acceptance Criteria

pH Value

Osmolality (1:10)

Incubation at 32°C

Incubation at 20°C

Accepted Result

7.3 - 7.5

n.a.

Sterile

Sterile

Composition

Inorganic Salts

Di-Sodium hydrogen phosphate	11,500.0000 mg/L
Potassium chloride	2,000.0000 mg/L
Potassium dihydrogen phosphate	2,000.0000 mg/L
Sodium chloride	80,000.0000 mg/L

Other Components

EDTA-Na ₂	2,000.0000 mg/L
Trypsin	5,000.0000 mg/L