

Specification

Product #: L 1861

Description: Trypsin powder (1 : 250) porcine origin, 25 g

Specification: Trypsin powder (1 : 250) porcine origin

Storage: 2°C - 8°C

Testings

Acceptance Criteria

Appearance

Solubility

Trypsin

Chymotrypsin

Protease Activity

Loss on drying

Accepted Result

white to whitish powder

partly soluble in water

≥225USP.U/mg

75 - 180USP.U/mg

225 - 275USP.U/mg

<5% vacuum oven at 60°C

Composition