

Specification

Product #: M 3104

Description: MEM NEAA (100x) w/o: L-Glutamine, 100 ml

Specification: w/o: L-Glutamine

Storage: 2°C - 8°C

Testings

Acceptance Criteria

pH Value

Osmolality

Incubation at 32°C

Incubation at 20°C

Accepted Result

1.3 - 1.5

n.a.

Sterile

Sterile

Composition

mg/L

Amino acids L-Alanine	890.0000 mg/L
Glycine	750.0000 mg/L
L-Asparagine x H ₂ O	1,500.0000 mg/L
L-Aspartic acid	1,330.0000 mg/L
L-Glutamic acid	1,470.0000 mg/L
L-Glutamine	0.0000 mg/L
L-Proline	1,150.0000 mg/L
L-Serine	1,050.0000 mg/L