

Specification

Product #: M 4601

Description: DPBS (10x) w/o: Ca and Mg, 500 ml

Specification: w/o: Ca w/o: Mg

Storage: 2° - Room Temperature

Testings

Acceptance Criteria

pH Value

Osmolality

Incubation at 32°C

Incubation at 20°C

Accepted Result

6.5 - 7.1

270 - 300 mOsm/kg

Sterile

Sterile

Composition

Inorganic Salts

di-Sodium hydrogen phosphate anhydrous	11,500.0000 mg/L
Potassium chloride	2,000.0000 mg/L
Potassium dihydrogen phosphate	2,000.0000 mg/L
Sodium chloride	80,000.0000 mg/L