

Specification

Product #: M 4611

Description: DPBS unsteril, w/o: Ca and Mg, 2,5 L

Specification: w/o: Ca w/o: Mg

Storage: 2° - Room Temperature

Testings

Acceptance Criteria

pH Value

Osmolality

Accepted Result

7.2 - 7.8

270 - 300 mOsm/kg

Composition

anhydrous

1150.00

1,150.0000 mg/L

Inorganic Salts

Potassium chloride

200.0000 mg/L

Potassium dihydrogen phosphate

200.0000 mg/L

Sodium chloride

8,000.0000 mg/L