

Specification

Product #: M 4621

Description: DPBS Powder 10L w/o: Ca and Mg, 10 L

Specification: DPBS Powder 10L w/o: Ca and Mg

Storage: 2° - Room Temperature

Testings

Acceptance Criteria

pH at RT (9.55 g/L)

Osmolality (9.55 g/L)

Accepted Result

7.2 - 7.8

262 - 290 mOsm/kg

Composition

Inorganic Salts

di-Sodium hydrogen phosphate anhydrous	1,150.0000 mg/L
Potassium chloride	200.0000 mg/L
Potassium dihydrogen phosphate	200.0000 mg/L
Sodium chloride	8,000.0000 mg/L