

Specification

Product #: M 4626

Description: DPBS w/o: Ca and Mg, 500 ml

Specification: w/o: Ca w/o: Mg

Storage: 2° - Room Temperature

Testings

Acceptance Criteria

pH Value

Osmolality

Incubation at 32°C

Incubation at 20°C

Accepted Result

7.2 - 7.8

270 - 300 mOsm/kg

Sterile

Sterile

Composition

Na₂HPO₄

1,150.0000 mg/L

Inorganic Salts

Potassium chloride

200.0000 mg/L

Potassium dihydrogen phosphate

200.0000 mg/L

Sodium chloride

8,000.0000 mg/L