

Specification

Product #: M 4906

Description: HBSS w/o: Ca and Mg, w: 0,35 g/l NaHCO₃, 500 ml

Specification: w/o: Ca and Mg w: 0.35 g/L NaHCO₃

Storage: 2° - Room Temperature

Testings

Acceptance Criteria

pH Value

Osmolality

Incubation at 32°C

Incubation at 20°C

Accepted Result

7.1 - 7.7

270 - 300 mOsm/kg

Sterile

Sterile

Composition

Inorganic Salts

Calcium chloride x 2H ₂ O	0.0000 mg/L
di-Sodium hydrogen phosphate dried	47.8800 mg/L
Magnesium sulphate anhydrous	0.0000 mg/L
Potassium chloride	400.0000 mg/L
Potassium dihydrogen phosphate	60.0000 mg/L
Sodium chloride	8,000.0000 mg/L

Other Components

D(+)-Glucose anhydrous	1,000.0000 mg/L
NaHCO ₃	350.0000 mg/L
Phenol red	10.0000 mg/L