

# #400024




## Cell Culture Plate 24 well

 award 2018

24 well Plate 85.4×127.5mm PS, sterile  
"plasma" growth area 1,9 cm<sup>2</sup> flat bottom  
Non - pyrogenic  
Non - cytotoxic  
DNase/RNase/DNA - free  
PS USP 6 SAL 10<sup>-6</sup>



Scan for COA



**LOT**  
201/968741

**EXP**  
6 Years

1 pc/bag

**STERILE** **R**

**Bioswisstec Ltd.** / e-mail : [info@bioswisstec.com](mailto:info@bioswisstec.com)  
website : [bioswisstec.com](http://bioswisstec.com) / Schaffhausen, Switzerland

**European Union Distribution / EC REP Laboglob.com GmbH**  
e-mail : [info@laboglob.com](mailto:info@laboglob.com) / website : [laboglob.com](http://laboglob.com)  
Gottmadingen, Germany